




Noninvasive & Continuous Hemoglobin (SpHb[®]) Monitoring

Lab Hb: Intermittent & Delayed Results While Transfusion Decisions Are Made in Real Time



How is Masimo's Solution Different?

Continuous SpHb provides real-time indication of *changes* in Hb, as well as when Hb is *stable*



How is Masimo's Solution Different?

With Masimo SpHb:


- **Noninvasive**
- **Continuous**



**Masimo SpHb permits
continuous measurement
of hemoglobin between
invasive blood samples**



rainbow® Multi-wavelength Technology



2 wavelengths of light used to measure oxygen saturation (SpO_2)

7+ wavelengths of light used by rainbow® to measure SpHb, SpCO, SpMet, and SpO_2



rainbow[®] Technology



SpO2, Pulse Rate plus *Additional Blood Constituents* which previously required a CO-Oximeter... is known as
“Pulse CO-Oximetry”




Multi-Parameter Trends of SpHb (and PVI) May Facilitate Transfusion Decisions




Intermittent Invasive Hemoglobin Sampling VS Real-time SpHb Trending Between Invasive Hb Samples

Intermittent Hb sampling may not provide the **full** picture...




- > Continuous SpHb provides real-time indication of whether hemoglobin is *stable*, *dropping*, or *rising*
- > Continuous visibility to Hb changes in between lab samples may help avoid unnecessary transfusions...
 - When *SpHb is stable* when Hb is perceived to be dropping, **or**
 - When *SpHb is rising* when Hb is perceived not to be rising fast enough, **AND...**
- > Detect post-op bleeding **IF** *SpHb is dropping* when Hb should be stable

Continuous SpHb: Stable Trend May Facilitate Delaying Transfusions at Lower Hb Levels



Continuous SpHb: Rising Trend May Facilitate Delaying Transfusion of Additional Units



Liver transplant– 69 y/o female

— SpHb (g/dL) ▲ tHb CO-Ox (g/dL) ● tHb Coulter Counter (g/dL)

Hemoglobin Measurement Variability

between Two of the Same Model Lab Device

- Two IDENTICAL models of five different types of laboratory CO-Oximeters
- Device A vs. Device B of each model analyzed 72 consecutive blood samples from 12 pts

Manufacturer Model	Radiometer ABL-735	Bayer Rapidlab 860	Nova STP CCX 1	Roche Omni S	Instrumentation 682 CO-Oximeter
Bias (g/dL)	0.0	-0.3	-0.8	-0.4	0.4
1st Standard deviation (g/dL)	0.1	0.2	0.3	0.8	1.2
2nd Standard deviation (g/dL)	0.2	0.4	0.6	1.6	2.4

Two calibrated lab devices of the same model type
1st standard deviation = Avg of 0.5 g/dL; highest 1.2 g/dL




Hemoglobin Measurement Variability Between Two Different Model Lab Devices

Two different lab devices

- pHOx™ CO-Oximeter (Nova) COULTER® Ac-T diff2™ (Beckman Coulter)

Consecutive blood samples measured on both devices




N=471 samples from 33 patients


**Two calibrated lab
devices analyzing
sequential blood
draws varied by up
to 2 g/dL**

SpHb, Dept CO-Oximeter, Hemocue *Point-to-Point Accuracy vs Lab CO-Oximeter*


- > 471 hemoglobin measurements from 62 Surgical ICU patients
- > 3 Hb methods vs. reference Hb (central laboratory hematology analyzer, Sysmex XT2000i)
 - SpHb, satellite CO-Oximeter (Siemens RapidPoint 405), point-of-care device (HemoCue 301)



1.0 g/dL
(A_{rms})




1.1 g/dL
(A_{rms})




1.3 g/dL
(A_{rms})

SpHb, Dept CO-Oximeter, Hemocue *TREND Accuracy vs Lab CO-Oximeter*


- > 471 hemoglobin measurements from 62 Surgical ICU patients
- > Changes in 3 Hb methods vs. changes in reference Hb (central laboratory hematology analyzer, Sysmex XT2000i)
 - SpHb, satellite CO-Oximeter (Siemens RapidPoint 405), point-of-care device (HemoCue 301)



R=0.64



R=0.60



R=0.39

**Continuous SpHb is shown to help reduce
transfusions – without changes in
transfusion protocol**



Randomized Controlled Trial in Lower Blood Loss Surgery

Objective, Patients, Randomization & Methods

> Objective

- Determine whether Continuous SpHb monitoring helps to reduce surgical transfusion frequency and average amount transfused in lower blood loss surgery

> Patients

- Orthopedic surgery (lower likelihood of transfusion) at academic medical center (Mass General Hospital, Boston, MA)

> Methods

- **Standard Care Group**


- Treat as normally would

- **SpHb Group**

- Treat as normally would but add continuous SpHb
 - **No transfusion protocol changes related to SpHb**




SpHb Monitoring Impact on Frequency of RBC Units Transfusions in Lower Blood Loss Surgery



Randomized controlled trial in 327 orthopedic surgery pts, 157-Standard Care & 170-SpHb

* p=0.03 vs. Standard Care Group;


SpHb Monitoring Impact on Average RBC Units Transfused per Patient in Lower Blood Loss Surgery



Randomized controlled trial in 327 orthopedic surgery pts, 157 Standard Care & 170 SpHb

* p=0.03 vs. Standard Care Group; † p=0.02 vs. Matched Retrospective Cohort

SpHb Monitoring Impact on Average RBC Units Transfused per Patient in Lower Blood Loss Surgery



Randomized controlled trial in 327 orthopedic surgery pts, 157 Standard Care & 170 SpHb

**p<0.0001 vs. Standard Care Group; †† p<0.0001 vs. Matched Retrospective Cohort

Randomized Controlled Trial

No Differences in Other Variables

- > **Frequency of patients receiving intraoperative Hb testing in SpHb and Standard Care Groups**
 - 11.8% vs. 16.3%, $p=ns$
- > **Mean number of Hb tests performed in SpHb and Standard Care Groups**
 - 0.24 vs. 0.21 tests per case, $p=ns$
- > **Safety variables**
 - No patient from either group received a transfusion during the immediate twelve-hour postoperative period
 - No differences at 28 days in the rate of post-operative complications between the SpHb and Standard Care Groups (1.9% vs. 3.0%, $p=ns$)

If SpHb Monitoring Reduces RBC Transfusions in Lower Blood Loss Surgery...

**Is it possible that SpHb monitoring
could have an even greater impact in
high blood loss surgery?**

Prospective Cohort Study in High Blood Loss Surgery

Objective, Patients, Randomization & Methods

> Objective

- Determine whether Continuous SpHb monitoring helps to reduce surgical transfusion frequency and average amount transfused in high blood loss surgery

> Patients

- Neurosurgery at academic medical center (Cairo University, Egypt)

> Methods

- **Standard Care Group**

- Treat as normally would


- **SpHb Group**

- Treat as normally would but add continuous SpHb

- **Both Groups**


- Blood samples taken at baseline and when EBL was $\geq 15\%$ of total blood volume
 - RBC transfusion initiated if hemoglobin was ≤ 10 g/dL and continued until the EBL was replaced and hemoglobin > 10 g/dL

SpHb Monitoring Impact on Frequency of ≥ 3 RBC Unit Transfusions in High Blood Loss Surgery



Prospective cohort study in 106 neurosurgery surgery pts, 61 Standard Care & 45 SpHb
* $p < 0.01$ vs. Standard Care Group

SpHb Monitoring Impact on Average RBC Units Transfused per Patient in High Blood Loss Surgery



Prospective cohort study in 106 neurosurgery surgery pts, 61 Standard Care & 45 SpHb
**p<0.001 vs. Standard Care Group



SpHb Clinical Benefits Extend Beyond Transfusion Decisions...

**Continuous SpHb can also help
identify post-op bleeding**




Detecting Post-Op Bleeding

- > **Bleeding is frequent in surgery, ICU, OB pts**
 - Up to 35% of patients¹
- > **Bleeding is a significant risk factor**
 - Late detection further increases the risk²
 - Responsible for 19% of in-hospital maternal deaths³
- > **Late bleeding detection impact on patient care**
 - Average hospital has multiple patients per year with serious injury or death due to late detection of bleeding⁴
- > **Bleeding significantly increases the total cost of patient treatment²**
- > **Low Hb identifies almost 90% of patients with bleeding⁵**
 - But traditional lab measurements are infrequent and delayed
- > **Joint Commission sentinel event alert for OB patients**
 - Calls for protocols to improve the ability to detect hemorrhage⁶



1 Hebert PC. Crit Care. 1999; 3(2):57-63. 2 Herwaldt LA. Infect Control Hosp Epidemiol. 2003; 24(1):44-50. 3 Bateman BT et al. Anesth Analg May 2010 110:1368-1373. 4 HRA Research of Hospital Executive 2012. 5 Bruns B et al. J Trauma. 2007; 63(2):312-5. 6 The Joint Commission, "Sentinel Event Alert: Preventing Maternal Death" Issue 44, January 26, 2010

Post-trauma Surgery Case Study



rainbow[®] Technology Clinical Benefits Extend Beyond Transfusion Decisions...

**PVI[®] helps with patient fluid management,
which impacts outcomes and
transfusion decisions**



Fluid Administration Challenges

- > **Fluid administration is critical**
 - To increase cardiac output and optimize patient status, enabling end organ preservation¹
- > **Unnecessary fluid administration associated with increased morbidity and mortality²**
- > **Fluid administration causes hemodilution**
 - Which reduces Hb concentration and increases likelihood of transfusion³
- > **Traditional static measurements not reliable to predict volume responsiveness**
 - CVP, SBP, DBP, Pulse Pressure, Wedge Pressure, Cardiac output⁴
- > **New dynamic monitoring technologies are effective**
 - Proven to improve outcomes but are invasive or complex - and costly⁵
 - Pulse pressure variation, stroke volume variation
 - Appropriate for higher-risk patients
- > **Opportunity to improve care in moderate to low risk population**
 - Using low cost, noninvasive technology



1 Perel A. Anesth Analg. 2008; 106 (4):1031-33 2 Bundgaard-Nielsen M et al. Acta Anaesthesiol Scand. 2007; 51(3):331-40

3 Huybregts RA et al. Anesth Analg. 2009. 109(2): 331-339. 4 Michard F et al. Chest. 2002; 121(6):2000-08 5 Joshi G et al. Anesth Analg. 2005; 101:601-5 .

Thank You

