

KHSC MHP Bedside Clinical Checklist for the TEAM CONTACT

To activate a Code Transfusion:

- Call 4444
- Provide the following information
 - Patient Location
 - Team contact's name
 - The patient's CR number
 - If the patient is in the OR or Connell 5 or going to the OR
 - Weight (for pediatric patients)
 - Biological Sex
 - Age
 - Type of Hemorrhage
 - History of antiplatelets or anticoagulants in the last 7 days



After the first blood pack arrives:

- Locate and carry the Code Transfusion phone (will be in the first pack)
- Ensure all necessary blood testing has been done (a pack of tubes will be sent with the first pack in a bag that has a red MHP sticker on it)
- Communicate directly with the Transfusion Medicine Laboratory (TML) (x4188) regarding additional blood products and lab results
- Direct the porter to the Core Lab to drop off samples or to the blood bank to pick up the next pack
- If there is a change in location or Team Contact notify the TML and hand over any remaining blood products and the Code Transfusion phone



Once MHP has been terminated:

- Inform the TML the MHP has been terminated
- Return all unused blood products in appropriate containers and Code Transfusion phone to the TML via the porter immediately**
- Inform the porter the MHP has been terminated so they can resume other tasks